




2022

**MJM**  
MEDICS

 Bseba, Baabda, Beirut, Lebanon

 [info@mjm-medics.com](mailto:info@mjm-medics.com)

 [www.mjm-medics.com](http://www.mjm-medics.com)

 +961 3168998 / +961 71660700







# Medical Supplies and Equipment

---

## Trusted Care

### About us

MJM group is a company founded in 2019, specialized in offering medical equipment and healthcare services. Our aim is to ultimately become the leading provider of the best medical services in Lebanon. We intend to do so by combining customer service and the implementation of the most advanced technologies in the medical field. We offer our customers superior technology and reliable quality of service.

As a distributor in medical health devices and equipment, we are constantly looking for new developments fulfilling the market needs. To support public and private clients in addressing the challenges of a complex and ever changing healthcare market, MJM Medics s.a.r.l. is able to create and implement innovative operational models which are sustainable, profitable, and whose vocation is to be adopted as standard worldwide.

MJM Medics s.a.r.l. is a young and dynamic company specializing in medical services, consumables and supplies.

We are offering our products based on: Reliability, Quality, Innovation and Service.

We are delivering innovative solutions for health care. Our products and concepts are developed in collaboration with medical professionals. These concepts must contribute in improving the quality of care for patients and professionals.

### Our Purposes

We are persistently to provide more developed health care, products, services, and assisting health care providers with their career, while adhering to our ethical values and respecting social and environmental concerns which will be positively reflected in our customers satisfaction, excellent teamwork performance and high-quality business sector.

### Our Missions

To be the biggest trading company in Lebanon in the medical field.

### Our Values

- Professionalism
- Creativity
- Trust
- Customer satisfaction
- Focusing on our client's advices
- Believe in our mission
- supporting leader to achieve vision

### Our Goals

- Projects commitment
- Respecting the timeframes
- Excellent financial outcomes
- Continues progress
- Patient safety
- Medical care provider development
- Fined the best solution for any medical problems



Medical Supplies and Equipment  
Trusted Care

## BASIC MEDICAL DEVICES

INFUSION THERAPY

UROLOGICAL CATHETERS AND URINE BAGS

OXYGEN THERAPY

SURGERY AND ANAESTHESIA

SYRINGES AND NEEDLES

BASIC DISPOSABLE EQUIPMENT

OTHER DISPOSABLE PRODUCTS

## DRESSINGS

TRADITIONAL DRESSING MATERIALS

ADHESIVE TAPES AND PLASTERS WITH DRESSING

SPECIAL PURPOSE DRESSINGS

BANDAGES

SPECIALIST PRODUCTS

## MEDICAL CLOTHING

COVERALL AND ISOLATION GOWN

SURGICAL GOWN

FACE MASKS

HEAD COVER AND SHOE COVER

# BASIC MEDICAL DEVICES



## INFUSION THERAPY

### I.V. cannula, sterile

#### Features

- type: pen type, with wing
- sizes: 14G, 16G, 18G, 20G, 22G, 24G, 26G
- high quality needle
- needle (stylet) made of stainless steel with back-cut blade
- handle that allows you to insert with one hand
- valve that prevents blood leakage to outside after inserting the cannula into the vessel
- flashback chamber secured with a luer lock cap
- latex-free
- phthalate-free
- EO sterilized



### Infusion Set

#### Features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with fluid containers/bags
- air inlet with anti-bacterial filter and coloured (blue) flap
- dropper: 20 drops = 1 ml +/- 0,1 ml
- fluid filter protecting against larger particles with pore size 15 µm
- 200 cm-long, soft, flexible tubing with luer lock conical fitting
- both tips of the device additionally secured with protective caps
- non-pyrogenic, non-toxic
- latex-free
- 'Y' injection port designed for additional injection.
- precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
- EO sterilized



### Infusion set with flow regulator

#### Features

- Air vented chamber with 15µm filter + roller clamp
- +needle free injection site + flow regulator + spin luer lock + prime stop cap
- designed to control flow rates from 5 ml/hr – 250 ml/hr manually.
- roller clamp for regulating fluid flow rate
- non-pyrogenic, non-toxic
- latex-free
- EO sterilized



### Extension tube for infusion pump

#### Features

- tubing made of PVC
- luer lock fitting - male and female
- latex-free, phthalate-free
- non-pyrogenic, non-toxic
- EO sterilized



## INFUSION THERAPY

### Transfusion set

#### Features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with bags
- precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
- air inlet with anti-bacterial filter and a coloured (red) flap
- dropper: 20 drops = 1 ml +/- 0,1 ml
- special blood filter with large area, pore size 200 µm
- 200 cm-long, soft, flexible tubing
- universal luer lock fitting at the end of tubing
- latex-free / non-pyrogenic, non-toxic
- EO sterilized



### 3-way stopcock

#### Features

- latex-free, phthalate-free / made of polycarbonate
- clear optical identifier of open/close position operating in 360 degrees
- consistent, straight flow channels and an optimal internal diameter ensure correct pressure measurement
- rotating male luer lock
- resistant to lipids
- non-pyrogenic, non-toxic
- EO sterilized



### 3-way stopcock with extension tube

#### Features

- clear optical identifier of open/close position
- three-way stopcock made of polycarbonate, extension tube made of PVC
- all inlets protected with caps
- resistant to lipids
- latex-free, phthalate-free
- non-pyrogenic, non-toxic
- EO sterilized



### Burette infusion set, Microdrip

#### Features

- 150 ml - Transparent, soft, cylindrical & calibrated measured volume chamber with bold marking scale.
- Automatic shut off valve prevents air trapping in the fluid line once the measured volume chamber is empty, and it also indicates the level of liquid.
- Sharp piercing spike provided with built in air inlet with bacterial barrier filter.
- Soft and kink resistant tubing.
- 'Y' injection port designed for additional injection.
- latex-free / non-pyrogenic, non-toxic
- EO sterilized



## UROLOGICAL CATHETERS AND URINE BAGS

### 2-way, silicone-coated latex, sterile

#### Features

- made of soft and flexible latex, that is resistant to bending and kinking
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- silicone-coated surface that allows atraumatic insertion of the catheter
- easy to fill and to empty, tear-resistant balloon of various capacities
- plastic or rubber valve to ensure balloon tightness
- contains latex, does not contain phthalates
- EO sterilized



### 2-way, 100% silicone, sterile

#### Features

- made of soft and flexible silicone, that reduces the allergic reaction of the mucous membranes and increases the patient's comfort
- atraumatic, soft rounded, special tip facilitating insertion of the catheter to the bladder
- tip with two lateral eyes with smooth edges
- radiopaque line along the entire length of the catheter
- easy to fill and to empty, tear-resistant balloon of various capacities
- plastic valve to ensure balloon tightness
- separately packed Catheter Spigot inside the packaging
- latex-free, phthalate-free
- EO sterilized



### Urine bag sterile

#### Features

- latex-free 2000 ml bag made of 0,22 mm thick PVC
- double sealed PVC for maximum protection against tearing the bag
- effective non-return valve
- tight T-valve, that can be easily operated with one hand
- precise, easy-to-read scale (every 100 ml), adequate for diuresis measurement
- tubing made of material that is resistant to bending and kinking and which ensures a free and effective urine outflow
- EO sterilized



### Urine bag with sample port 2W, sterile

#### Features

- latex-free 2000 ml bag made of medical grade PVC
- white back wall making urine easy to notice
- effective non-return valve
- tag to describe the patient's details and date
- additional fixing with a string
- self-tightening port for sample taking at the end of tubing
- Pasteur type droplet chamber providing a „dry“ air barrier
- at the place where the tubing connects to the chamber, additional reinforcing element prevents bending and kinking the tubing
- precise, easy-to-read scale, adequate for diuresis measurement, with high measurement accuracy
- tight T-valve, that can be easily operated with one hand and which can be folded and tucked in a hard sheath
- two hydrophobic air filters (in a bag and chamber), which prevent sucking, balancing the pressure inside the system
- 120 cm long, large diameter tubing, ensuring a free and effective urinary outflow with a universal stepped connector
- tubing equipped with a sliding clamp and a fastening pin to stabilize the tubing
- integrated universal hanger system matching the round and square bed rails
- EO sterilized



## OXYGEN THERAPY

### Oxygen mask

#### Features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- equipped with a 210 cm long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- revolving connector that allows adjusting to patient's position
- latex-free
- EO sterilized



### Non-rebreather mask

#### Features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- equipped with a 210 cm long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- 1000 ml oxygen reservoir
- revolving connector that allows adjusting to patient's position
- silicone valves on the connector and ventilation openings
- latex-free
- EO sterilized



### Nebulizer mask

#### Features

- made of transparent, non-toxic PVC
- adjustable nose metal plate and rubber fastening
- equipped with a 210 cm long tubing with universal connectors
- tubing with star cross-section that is resistant to bending
- 8 ml nebulizer (graduated every 1 ml)
- average size of the molecule at a flow of 6l / min - 3,23  $\mu$ m
- latex-free
- EO sterilized



### Nasal oxygen cannula

#### Features

- made of medical PVC
- very soft tips with smooth edges
- fastening under the chin with a soft, sliding adjustment ring
- various spacing of the nasal part of the cannula depending on the size
- universal connector compatible with any source of oxygen
- tubing with star cross-section that is resistant to bending
- latex-free
- EO sterilized



## OXYGEN THERAPY

### Oxygen connecting tube

#### Features

- latex-free / made of transparent, non-toxic PVC
- inner ribbed design to avoid choking
- tubing with star cross-section that is resistant to bending, with reinforcing straps along the entire length
- standardized connectors at the end of the tubing allowing connection to oxygen devices
- EO sterilized



### Lung Exerciser - 3 Ball Spirometer

#### Features

- color-coded balls in each chamber prove to be visual incentive to the patient.
- connecting tube with 12mm OD connector and mouth piece.
- flow rates 600ml/sec, 900ml/sec, and 1200ml/sec by using different colors of ball for easy identification of the flow rates.
- can be dis-assembled into parts for cleaning and disinfection.
- non-sterile, individually packed in a box.
- compact, high quality break-resistant plastic, 3-ball incentive spirometer.



## SURGERY AND ANAESTHESIA

### Endotracheal tube uncuffed, sterile

#### Features

- latex-free made of flexible PVC
- radiopaque line over the entire tube length
- equipped with lateral Murphy eye
- graduated every 1 cm
- double depth marker to control the position of the tube
- tube size marked in 3 places (on the connector and in two places on the tube body)
- detachable connector with an outer diameter of 15 mm
- EO sterilized



### Endotracheal tube cuffed, sterile

#### Features

- latex-free made of flexible PVC
- equipped with lateral Murphy eye and radiopaque line over the entire tube length
- graduated every 1 cm, double depth marker to control the position of the tube
- low-pressure, high-volume cuff
- tube size marked in 4 places (on the connector, on the pilot balloon and in two places on the tube body)
- detachable connector with an outer diameter of 15 mm
- pilot balloon with LOT and size indication
- EO sterilized



## SURGERY and ANAESTHESIA

### Carbon Steel Surgical

#### Features

- sterile surgical carbon(stainless) steel blades.
- used for the basic surgical surgery for cutting soft tissues
- prevent the cross infections from the medical sources and suitable to the field rescue
- way of sterilization sterilized by gamma radiation.



## SYRINGES and NEEDLES

### Disposable hypodermic needle

#### Features

- needle tube is made of smooth and durable stainless steel
- sheath and hub made of polypropylene
- thin walls assure adequate stream pace
- tip sharpened in three surfaces (back cut)
- bevel type: LB/BL standard long bevel
- surface of needle tube covered with sliding agent (silicone)
- glued needle elements to prevent their separation (glue - epoxy resin)
- hubs coloured according to ISO code allowing fast recognition of needle size
- latex-free, phthalate-free / non-toxic, non-pyrogenic
- EO sterilized



### Three Parts Disposable Syringe With and Without Needle Luer Lock Tip

#### Features

- material: medical grade high transparent PP
- central nozzle, Luer lock tip
- product structure: barrel, plunger, piston and hypodermic needle
- individual packing: PE bag, blister package
- inner packing: inner box (card paper or corrugated paper)
- outer packing: corrugated carton
- sterilized by EO gas / non-toxic, non-pyrogenic



### Insulin Syringe

#### Features

- with bypacked needle 0,4 x 13 mm
- transparent barrel made of polypropylene allows full visua control of content
- sealing in the form of double ring on plunger stopper situated on the top of plunger
- smooth plunger sliding thanks to rubber sealing
- limitation flange securing against accidental removal of plunger
- special design of plunger stopper reduces the residual volume of medicine/liquid
- clear, legible and permanent black graduation facilitating dosage
- latex-free, phthalate-free / non-pyrogenic, non-toxic
- EO sterilized



## BASIC DISPOSABLE EQUIPMENT

### Suction catheter

#### Features

- made of soft and flexible PVC, that is resistant to bending and kinking
- frozen surface enabling easy insertion
- atraumatic, soft rounded, open tip
- tip with two opposite lateral eyes with smooth edges
- colour-coded connectors for size recognition, numerical size
- latex-free
- phthalate-free
- EO sterilized



### Yankauer surgical suction cannula

#### Features

- latex-free
- made of soft and flexible PVC
- along the tube a special longitudinal reinforcement to prevent bending and sucking of the tube
- perfectly smooth inner surface of tube prevents the sucked content from settling on the walls, ensuring correct and optimal flow
- connector type funnel-funnel cut-to-fit with the possibility of cutting in the marked places of the connector and matching to any type of medical suction device
- funnel connector equipped with a special bending system, enabling convenient manipulation of the tube
- tube connectors compatible with the standard suction cannulas, enabling their reliable, effective and tight connection
- EO sterilized



### vaginal speculum

#### Features

- both blades made of transparent polystyrene
- locking bolt made of polyethylene
- blades position blocked by a quarter turn of the locking bolt (without screwing in the bolt)
- perfectly smooth and streamlined edges
- colour coded bolts depending on the size
- single-use
- latex-free, phthalate-free
- EO sterilized



## BASIC DISPOSABLE EQUIPMENT

### Wound Drainage Reservoir

#### Features

- 3-spring Evacuator provides reliable suction for surgical cases with larger amounts of drainages.
- transparent Sidewalls and Clear Volume Calibrations facilitate identification and measurement of exudate.
- drain Port with Attached Reclosable Pouring Spout is provided permits convenient emptying of the reservoir unit.
- suction Port Equipped with an Anti-reflux Valve to preclude retrograde flow of drainage fluid during emptying and reactivation.
- available in convenient kit or separately with round or flat drains, with pre-attached trocar or not for procedures flexibility.
- sterile, for single-use only
- superior biocompatibility, flexibility and great comfort to patient, which delivers extremely low tissue response to patient, which deliver extremely low tissue; response to wound, which minimize possible tissue adhesion.
- universal stepped adapter allows connecting to kinds of suction tube.
- embedded X-ray thread enables easy detection and position of suction tube.
- autotraumatic tip less pain of patient.
- anti-reflux valve totally eliminates liquid reflux.



## OTHER DISPOSABLE PRODUCTS

### Wooden tongue depressor

#### Features

- made of birch wood
- rounded edges preventing injuries
- natural colour
- non-sterile



### PE gloves

#### Features

- made of polyethylene
- coarse surface and tear-resistant weld
- ambidextrous
- elongation min. 130%, tensile strength min. 10 MPa
- approved for contact with food



### Clinical Thermometer Mercury

#### Features

- medical Clinical Mercury Thermometer For Rectal Use
- clinical thermometer is made of glass and mercury
- measurement range: 35-42 in C scale or 94-108 in F scale
- available usage: rectal
- accuracy: -0.15C ~ +0.1C



# DRESSINGS



## TRADITIONAL DRESSINGS MATERIALS

### Non sterile Gauze Swabs

#### Features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- edges ES-type folded, which prevents fraying of loose threads
- available with a sewed in X-ray detectable thread, which enables fast and non-invasive identification of the swab left in the operating field



without X-ray detectable thread

dimensions	no. of threads	no. of plies	packaging unit paper bag
5 cm x 5 cm	13	8	100 pcs
7,5 cm x 7,5 cm	13	8	100 pcs
10 cm x 10 cm	13	8	100 pcs
5 cm x 5 cm	13	12	100 pcs
5 cm x 5 cm	13	12	500 pcs
7,5 cm x 7,5 cm	13	12	100 pcs
7,5 cm x 7,5 cm	13	12	500 pcs
10 cm x 10 cm	13	12	100 pcs
10 cm x 10 cm	13	12	500 pcs
5 cm x 5 cm	17	8	100 pcs
7,5 cm x 7,5 cm	17	8	100 pcs
10 cm x 10 cm	17	8	100 pcs
10 cm x 20 cm	17	8	100 pcs
5 cm x 5 cm	17	12	100 pcs
5 cm x 5 cm	17	12	500 pcs
7,5 cm x 7,5 cm	17	12	100 pcs
7,5 cm x 7,5 cm	17	12	500 pcs
10 cm x 10 cm	17	12	100 pcs
10 cm x 10 cm	17	12	500 pcs
10 cm x 20 cm	17	12	100 pcs
5 cm x 5 cm	17	16	100 pcs
7,5 cm x 7,5 cm	17	16	100 pcs
10 cm x 10 cm	17	16	100 pcs
10 cm x 20 cm	17	16	100 pcs

with X-ray detectable thread

dimensions	no. of threads	no. of plies	packaging unit paper bag
5 cm x 5 cm	17	8	100 pcs
7,5 cm x 7,5 cm	17	8	100 pcs
10 cm x 10 cm	17	8	100 pcs
5 cm x 5 cm	17	12	100 pcs
7,5 cm x 7,5 cm	17	12	100 pcs
10 cm x 10 cm	17	12	100 pcs
10 cm x 20 cm	17	12	100 pcs
5 cm x 5 cm	17	16	100 pcs
7,5 cm x 7,5 cm	17	16	100 pcs
10 cm x 10 cm	17	16	100 pcs
10 cm x 20 cm	17	16	100 pcs

### Sterile Lap Sponge

#### Features

- 100% hydrophilic cotton
- bleached with the use of hydrogen peroxide
- with a sewed in X-ray detectable thread or chip, which enables fast and non-invasive identification of the lap sponge left in the operating field
- with a belt to fix the lap sponge beyond the operating field
- Double wapping
- Pre-washed
- 45 cm x 45 cm
- 13 threads
- with blue loop



## ADHESIVE TAPES AND PLASTERS WITH DRESSING

### Features

- flexible, conforms to the body's contours, guarantees freedom of movement
- characterized by adequate air and water vapour permeability
- good adhesion
- can be divided lengthwise and widthwise without the use of scissors
- hypoallergenic

### Surgical Fabric Tape

It is used to fix medical devices on skin.

100% Fabric, Zinc oxide adhesive, Non sterile

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V1101	9.14 m x 1.25 cm	24 Pcs/Box	10 boxes / 240 pcs
V1102	9.14 m x 2.5 cm	12 Pcs/Box	10 boxes / 120 pcs
V1105	9.14 m x 5 cm	6 Pcs/Box	10 boxes / 60 pcs
V1107	9.14 m x 7.5 cm	4 Pcs/Box	10 boxes / 40 pcs
V1110	9.14 m x 10 cm	2 Pcs/Box	10 boxes / 20 pcs



### Surgical Silk Tape

It is used to fix medical devices on skin.

Silk-like fabric, Zinc oxide adhesive, Non sterile

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V2101	9.14 m x 1.25 cm	24 Pcs/Box	10 boxes / 240 pcs
V2102	9.14 m x 2.5 cm	12 Pcs/Box	10 boxes / 120 pcs
V2105	9.14 m x 5 cm	6 Pcs/Box	10 boxes / 60 pcs
V2107	9.14 m x 7.5 cm	4 Pcs/Box	10 boxes / 40 pcs
V2110	9.14 m x 10 cm	2 Pcs/Box	10 boxes / 20 pcs



### Surgical Transparent Tape

It is used to fix medical devices on skin. It is suitable for use in treatment areas which should be monitored visually.

Pe, Perforated, Hypoallergenic, Breathable, Non-sterile

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V3101	9.14 m x 1.25 cm	24 Pcs/Box	10 Boxes / 240 Pcs
V3102	9.14 m x 2.5 cm	12 Pcs/Box	10 Boxes / 120 Pcs
V3105	9.14 m x 5 cm	6 Pcs/Box	10 Boxes / 60 Pcs
V3107	9.14 m x 7.5 cm	4 Pcs/Box	10 Boxes / 40 Pcs
V3110	9.14 m x 10 cm	2 Pcs/Box	10 Boxes / 20 Pcs



### Surgical Paper Tape

It is used to fix medical devices on skin. Especially suitable for sensitive skin.

Paper, Hypoallergenic, Conformable and highly breathable, Non-sterile

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V5101	9.14 m x 1.25 cm	24 Pcs/Box	10 Boxes / 240 Pcs
V5102	9.14 m x 2.5 cm	12 Pcs/Box	10 Boxes / 120 Pcs
V5105	9.14 m x 5 cm	6 Pcs/Box	10 Boxes / 60 Pcs
V5107	9.14 m x 7.5 cm	4 Pcs/Box	10 Boxes / 40 Pcs



## SPECIAL PURPOSE DRESSINGS

### Sterile NON WOVEN Wound Dressing

Suitable for all kind of skin injuries and treatments. It's used for protecting post operative wound area and sutured skin. The person can make his/her own dressing easily. Sterile, Non Woven, Hypoallergenic, Breathable

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V607	5 cm x 7.5 cm	50 Pcs/Box	4 Boxes / 200 Pcs
V609	5 cm x 9 cm	50 Pcs/Box	4 Boxes / 200 Pcs
V610	9 cm x 10 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V615	9 cm x 15 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V620	9 cm x 20 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V625	9 cm x 25 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V630	9 cm x 30 cm	25 Pcs/Box	4 Boxes / 100 Pcs



### Sterile Transparent PU Wound Dressing

Suitable for all kind of skin injuries and treatments. It's used for protecting post operative wound area and sutured skin. The person can make his/her own dressing easily. Sterile, Pu, Hypoallergenic, Breathable

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V6107P	5 cm x 7.5 cm	50 Pcs/Box	4 Boxes / 200 Pcs
V6109P	5 cm x 9 cm	50 Pcs/Box	4 Boxes / 200 Pcs
V6110P	9 cm x 10 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V6115P	9 cm x 15 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V6120P	9 cm x 20 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V6125P	9 cm x 25 cm	25 Pcs/Box	4 Boxes / 100 Pcs
V6130P	9 cm x 30 cm	25 Pcs/Box	4 Boxes / 100 Pcs



### Sterile NON WOVEN I.V. Fixation Dressing

It is used to fix the catheter to the body. Sterile ,Non Woven , Hypoallergenic, Breathable.

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V6608	6 cm x 8 cm	50 Pcs/Box	20 Boxes / 1000 Pcs
V6608P	6 cm x 8 cm	50 Pcs/Box	20 Boxes / 1000 Pcs
V6710	7.2 cm x 10 cm	50 Pcs/Box	20 Boxes / 1000 Pcs



### Sterile PU Transparent I.V. Fixation Dressing

It is used to fix the catheter to the body. Sterile , Transparent Pu, Hypoallergenic, Breathable

Ref No	Dimentions	Pcs/Package	Pcs/Carton
V6608F	6 cm x 8 cm	50 Pcs/Box	20 Boxes / 1000 Pcs
V6608FP	6 cm x 8 cm	50 Pcs/Box	20 Boxes / 1000 Pcs
V6710F	7.2 cm x 10 cm	50 Pcs/Box	20 Boxes / 1000 Pcs



## SPECIAL PURPOSE DRESSINGS

### Surgical Elastic Fixation Tape

It is used for fixing the medical devices in flexible and active areas on the skin.  
Non Woven, Hypoallergenic, Breathable, Non-sterile

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V4102	10 m x 2.5 cm	1 Pcs/Box	200 individual boxes
V4105	10 m x 5 cm	1 Pcs/Box	120 individual boxes
V4110	10 m x 10 cm	1 Pcs/Box	60 individual boxes
V4115	10 m x 15 cm	1 Pcs/Box	40 individual boxes
V4120	10 m x 20 cm	1 Pcs/Box	32 individual boxes
V4130	10 m x 30 cm	1 Pcs/Box	20 individual boxes



### Standard First Aid Plaster

Pe, water resistant, breathable.

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V855	19 mm x 72 mm	100 Pcs/Box	10 Boxes



### Injection Plaster

It is used to close the area after all kind of injection process.  
Pe, Perforated, Breathable, Hypoallergenic, Skin Color

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V8025	25 mm	100 Pcs/Box	100 Boxes
V8125	25 mm	1000 Pcs/Box	10 Boxes



### Sterile Hypoallergenic Eye Pad

It can be used for the treatment of injuries, allergic reactions and post-operative wounds on eye.  
Sterile, White Non Woven, Hypoallergenic, Breathable

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V6505	5 cm x 6.2 cm	25 Pcs/Box	12 Boxes / 300 Pcs
V6825	5.7 cm x 8.2 cm	25 Pcs/Box	12 Boxes / 300 Pcs
V6950	6.5 cm x 9.5 cm	50 Pcs/Box	10 Boxes / 500 Pcs



## BANDAGES

### Elastic Bandage

200% Elastic, Breathable, Washable and reusable, Anti Allergenic

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V80645	6 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	600 Pcs
V80845	8 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	500Pcs
V81045	10 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	400 Pcs
V81245	12 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	300 Pcs
V81545	15 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	250 Pcs
V82045	20 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	200 Pcs



### Cotton Crepe Bandage

100% Elastic, Breathable, Washable and reusable, Anti Allergenic.

RefNo	Dimintions	Pcs/Package	Pcs/Carlton
V20645	6 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	600 Pcs
V20845	8 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	500 Pcs
V21045	10 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	400 Pcs
V21245	12 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	300 Pcs
V21545	15 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	250 Pcs
V22045	20 cm x 4,5 m (stretched)	10 Pcs/Pe Pack	200 Pcs



### Gauze Bandage

100% Cotton, Breathable, Hydrophilic, Non Sterile

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V60540	5 cm x 4 m	50 Pcs/Pe Pack	2250 Pcs
V60740	7 cm x 4 m	50 Pcs/Pe Pack	1500 Pcs
V61040	10 cm x 4 m	50 Pcs/Pe Pack	1200 Pcs
V61540	15 cm x 4 m	25 Pcs/Pe Pack	750 Pcs
V62040	20 cm x 4 m	25 Pcs/Pe Pack	650 Pcs



### Elastic Gauze Bandage

100% Elastic, Breathable, Soft & Light, Non Sterile

Ref No	Dimintions	Pcs/Package	Pcs/Carlton
V70540	5 cm x 4 m (stretched)	50 Pcs/Pe Pack	2250 Pcs
V70740	7 cm x 4 m (stretched)	50 Pcs/Pe Pack	1500 Pcs
V71040	10 cm x 4 m (stretched)	50 Pcs/Pe Pack	1200 Pcs
V71540	15 cm x 4 m (stretched)	25 Pcs/Pe Pack	750 Pcs
V72040	20 cm x 4 m (stretched)	25 Pcs/Pe Pack	650 Pcs




## SPECIALIST PRODUCTS

### Paraffin Gauze Dressing

#### Features

- wide-mesh dressing made of cotton gauze covered with paraffin
- structure of the carrier material allows exudate to be freely drained to the covering secondary dressing
- covering with a layer of paraffin reduces risk of adhesion to the wound and allows a painless change of dressing
- neutral properties of paraffin allow for using the dressing in combination with other preparations containing medicinal substances or antiseptics
- can be shaped and cut according to the needs aseptically
- non-adhesive, necessary fixing with a tape, bandage or net tubular dressing
- hypoallergenic



### Alcohol Swab

#### Features

- made of non-woven
- soaked with isopropyl alcohol (70%)
- two pouches connected in one-piece packaging can be separated thanks to perforation



### Pharmaplast Incifilm Drape - Iodine Theatre Incise Drape.

#### Features

- frame border gives support to the edges of Difa Iodine Frame style, preventing wrinkles during application.
- adheres right up to the wound edge.
- yellowish Transparent surface.
- thin but extremely strong.
- conformable and elastic.
- provides a sterile working field around the operation site and a broad spectrum antimicrobial activity.
- holds the drapes securely in place.
- sterile and individually wrapped in peel apart pouch.
- sterilized by gamma irradiation.
- bacterial proof and waterproof film.
- available in a variety of sizes to meet your procedural and incision-size needs.



# MEDICAL CLOTHING




## COVERALL and ISOLATION GOWN

### Protective Coverall

#### Features

- made of a laminate (non-woven polypropylene and polyethylene) with a weight of 70 g/m<sup>2</sup>
- polyester seams additionally covered with tape
- equipped with a hood with elastic trim
- equipped with an elastic waist band, wrists and ankles
- zipper covered with a self-adhesive flap
- available in version with boots

SIZE ROZMIAR						
cm	86-94	94-102	102-110	110-118	118-129	129-141
158-166	S					
166-174		M				
174-182			L			
182-190				XL		
190-198					XXL	
198-206						XXXL

### Isolation Gown with elastic cuffs

#### Features

- made of polypropylene non-woven 20 g/m<sup>2</sup>, 25 g/m<sup>2</sup>, 30 g/m<sup>2</sup> or 40 g/m<sup>2</sup>
- ties at waist and neck
- long sleeves with elastic cuffs
- breathable
- available sizes: L, XL, XXL
- single-use



### Isolation Gown with knitted cuffs

#### Features

- made of polypropylene non-woven 20 g/m<sup>2</sup>, 25 g/m<sup>2</sup>, 30 g/m<sup>2</sup> or 40 g/m<sup>2</sup>
- ties at waist and neck
- long sleeves with knitted cuffs
- breathable
- available sizes: L, XL, XXL
- single-use



## SURGICAL GOWN

### Sterile Surgical Gown Standard





size	trimming	packaging unit
M		1 pc.
L		1 pc.
XL		1 pc.
XXL		1 pc.

includes  
2 absorbent towels  
30cm x 30cm

45 g/m<sup>2</sup>



### Sterile Surgical Gown Reinforced

size	trimming	packaging unit
M		1 pc.
L		1 pc.
XL		1 pc.
XXL		1 pc.

includes  
2 absorbent towels  
30cm x 30cm

45 g/m<sup>2</sup>



## STANDARDS GUIDE

CHARACTERISTICS	STANDARD	UNIT	REQUIREMENT			
			STANDARD REQUIREMENT		HIGH REQUIREMENT	
			CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA	CRITICAL PRODUCT AREA	LESS CRITICAL PRODUCT AREA
microbial penetration - dry	EN ISO 22612	CFU	not required	≤ 300 <sup>a</sup>	not required	≤ 300 <sup>a</sup>
microbial penetration - wet	EN ISO 22610	I <sub>b</sub>	≥ 2,8 <sup>b</sup>	not required	6,0 <sup>bc</sup>	not required
cleanliness - microbial bioburden	EN ISO 11737-1	CFU/100 cm <sup>2</sup>	≤ 300	≤ 300	≤ 300	≤ 300
particulate release	EN ISO 9073-10	IPM	≤ 4,0	≤ 4,0	≤ 4,0	≤ 4,0
linting	EN ISO 9073-10	log <sub>10</sub> (lint count)	≥ 40	≥ 40	≥ 40	≥ 40
liquid penetration	EN ISO 811	cm H <sub>2</sub> O	≥ 30	≥ 10	≥ 100	≥ 10
bursting strength - dry	EN ISO 13938-1	kPa	≥ 40	≥ 40	≥ 40	≥ 40
bursting strength - wet	EN ISO 13938-1	kPa	≥ 40	not required	≥ 40	not required
tensile strength - dry	EN 29073-3	N	≥ 20	≥ 20	≥ 20	≥ 20
tensile strength - wet	EN 29073-3	N	≥ 20	not required	≥ 20	not required

## FACE MASKS

### Face Mask

#### Features

- made of three layers of dust-free and hypoallergenic non-woven
- modelling nose-piece
- splash resistant
- protects operative field against exfoliated epidermis, hair and droplet infections
- provides effective protection in patients against infections transmitted from personnel's mouths and noses
- packaging: ergonomic cardboard box to remove individual pieces
- single-use



With Ties



With Earloops



With Headloops

### KN95 Mask

#### Features

- 3D Foldable Design
- Particle Filtering Efficiency (PFE)  $\geq 95\%$
- Elastic Ear-loop
- Adjustable and Invisible Nose Clip
- 5-Layer Protection (high quality nonwoven fabric, hot air cotton and melt-blown fabric)
- Soft Fabric and Comfortable to Wear
- single-use



## HEAD COVER and SHOE COVER

### Medical Head Cover

#### Features

- cap with an elastic rubber band
- diameter: approx. 53 cm
- made of airy polypropylene non-woven
- excellent breathability and wearing comfort
- packaging: ergonomic cardboard box to remove individual pieces
- single-use, non-sterile



### Shoe Cover

#### Features

- made of polypropylene non-woven 40 g/m<sup>2</sup>
- for shoe covers with anti-skid layer)
- elastic rubber band at the opening
- protection against contaminating rooms where higher cleanliness standards are required
- increased durability
- single-use





Medical Supplies and Equipment

---

*Trusted Care*



Bseba, Baabda, Beirut, Lebanon



[info@mjm-medics.com](mailto:info@mjm-medics.com)



[www.mjm-medics.com](http://www.mjm-medics.com)



+961 3168998 / +961 71660700