

## XNY GUIDANCEPORT™ DISPOSABLE TROCAR

Guidanceport™  
ENDOSCOPIC ACCESS

XNY's company moto of "CONSTANT INNOVATION" means customers receive exceptional performing devices which offer unique benefits that are functional and ergonomic. The drive to improve our products means greater levels of safety for patients and positive user experiences.



**XNY MEDICAL**  
Smart Surgical, Smart People, Smart Life

**NOVA**  
SURGICAL STAPLERS

**Guidanceport™**  
ENDOSCOPIC ACCESS

Address: N.51 Shuishan Rd., Zhonglou Economic Development District, Changzhou, Jiangsu, China, Zip Code:213023

TEL: +86-519-83272327

Fax: +86-519-86865852

E-mail: info@xnystaplers.com

Website: www.xnymedical.com





Guidanceport™ Disposable Bladeless Trocar



Guidanceport™ Disposable Bladed Trocar




Guidanceport™ Disposable Spiral Trocar




Guidanceport™ Disposable Optical Trocar

**ROBUST AND STRONG  
UNIVERSAL SEAL**  
Accommodates 5-12mm  
instruments without  
leakage



**UNIQUE UNIVERSAL  
SEAL DESIGN**  
Prevent seal damage  
and leakage of gas



**FUNNEL LEAD IN**  
Assists easy guided  
instrument insertion




**SPIRAL TIP**  
Less damage to  
abdominal wall  
and other tissues


**TRADITIONAL  
TISSUE  
SEPARATING TIP**  
Sizes range from  
3mm to 15mm,  
enables entry of  
instruments of  
different diameters



**TRANSPARENT TIP**  
Controllable insertion  
under endoscope



**DILATING TIP**  
Fast retraction of blade  
upon insertion that  
prevents further damage



## XNY GUIDANCEPORT™ DISPOSABLE TROCAR

**Guidanceport™**  
ENDOSCOPIC ACCESS

Guidanceport™ Disposable Bladeless Trocar	Code	Size	Colour	Pack Qty
	D3×100	φ4mm*100mm	White	6sets/box
	D5×100	φ6mm*100mm	Blue	6sets/box
	D10×100	φ11mm*100mm	Yellow	6sets/box
	D12×100	φ13mm*100mm	Green	6sets/box
	D15×100	φ16mm*100mm	Black	6sets/box
	D5×150	φ6mm*150mm	Blue	6sets/box
	D10×150	φ11mm*150mm	Yellow	6sets/box
	D12×150	φ13mm*150mm	Green	6sets/box
	D3A×100	φ4mm*100mm	White	6sets/box
	D5A×100	φ6mm*100mm	Blue	6sets/box
	D10A×100	φ11mm*100mm	Yellow	6sets/box
	D12A×100	φ13mm*100mm	Green	6sets/box
	D15A×100	φ16mm*100mm	Black	6sets/box
	D5A×150	φ6mm*150mm	Blue	6sets/box
	D10A×150	φ11mm*150mm	Yellow	6sets/box
	D12A×150	φ13mm*150mm	Green	6sets/box

Guidanceport™ Disposable Optical Trocar	Code	Size	Colour	Pack Qty
	DO5×100	φ6mm*100mm	Blue	6sets/box
	DO10×100	φ11mm*100mm	Yellow	6sets/box
	DO12×100	φ13mm*100mm	Green	6sets/box
	DO5×150	φ6mm*150mm	Blue	6sets/box
	DO10×150	φ11mm*150mm	Yellow	6sets/box
	DO12×150	φ13mm*150mm	Green	6sets/box
	DO5A×100	φ6mm*100mm	Blue	6sets/box
	DO10A×100	φ11mm*100mm	Yellow	6sets/box
	DO12A×100	φ13mm*100mm	Green	6sets/box
	DO5A×150	φ6mm*150mm	Blue	6sets/box
	DO10A×150	φ11mm*150mm	Yellow	6sets/box
	DO12A×150	φ13mm*150mm	Green	6sets/box

Guidanceport™ Disposable Spiral Trocar	Code	Size	Colour	Pack Qty
	DS5×100	φ6mm*100mm	Blue	6sets/box
	DS10×100	φ11mm*100mm	Yellow	6sets/box
	DS12×100	φ13mm*100mm	Green	6sets/box
	DS5A×100	φ6mm*100mm	Blue	6sets/box
	DS10A×100	φ11mm*100mm	Yellow	6sets/box
	DS12A×100	φ13mm*100mm	Green	6sets/box

Guidanceport™ Disposable Bladed Trocar	Code	Size	Colour	Pack Qty
	DK5×100	φ6mm*100mm	Blue	6sets/box
	DK10×100	φ11mm*100mm	Yellow	6sets/box
	DK12×100	φ13mm*100mm	Green	6sets/box
	DK5A×100	φ6mm*100mm	Blue	6sets/box
	DK10A×100	φ11mm*100mm	Yellow	6sets/box
	DK12A×100	φ13mm*100mm	Green	6sets/box

Guidanceport™ Disposable Trocar Cannula	Code	Size	Colour	Pack Qty
	E3	φ4mm*100mm	White	6pcs/box
	E5	φ6mm*100mm	Blue	6pcs/box
	E10	φ11mm*100mm	Yellow	6pcs/box
	E12	φ13mm*100mm	Green	6pcs/box
	E15	φ16mm*100mm	Black	6pcs/box

Veress Insufflation Needle	Code	Size	Colour	Pack Qty
	VN120	φ2mm*120mm		10 pcs/box
	VN150	φ2mm*150mm		10 pcs/box